

BH-Vis

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BH-Vis is a dual-photon fluorescent probe that targets the cell plasma membrane (CPM). It demonstrates significant potential for accurately identifying the levels of cellular apoptosis during the development of abdominal aortic aneurysm (AAA) in a mouse model.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481