

Dxd-D5

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dxd-D5 isomerase I inhibitor with an IC50 value of 0.31 $\mu$ M. It serves as the effective payload for the antibody-drug conjugate (ADC) targeting HER2, known as DS-8201a.
Targets(IC50)	ADC Cytotoxin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481