


Bongardol

Chemical Properties

CAS No. :	123690-76-6
Formula:	C36H64O3
Molecular Weight:	544.89
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bongardol is a natural product for research related to life sciences. The catalog number is TN3529 and the CAS number is 123690-76-6.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8352 mL	9.1762 mL	18.3523 mL
5 mM	0.367 mL	1.8352 mL	3.6705 mL
10 mM	0.1835 mL	0.9176 mL	1.8352 mL
50 mM	0.0367 mL	0.1835 mL	0.367 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481