

2E-Decenoylpiperidide

Chemical Properties


CAS No. :

Formula: C₁₅H₂₇NO

Molecular Weight: 237.387

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2E-Decenoylpiperidide is a useful organic compound for research related to life sciences and the catalog number is T125821.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2125 mL	21.0624 mL	42.1248 mL
5 mM	0.8425 mL	4.2125 mL	8.425 mL
10 mM	0.4212 mL	2.1062 mL	4.2125 mL
50 mM	0.0842 mL	0.4212 mL	0.8425 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481