


m-PEG4-sulfonic acid

Chemical Properties

CAS No. :	787524-78-1
Formula:	C9H20O7S
Molecular Weight:	272.32
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	m-PEG4-sulfonic acid is a PEG-based linker for PROTACs that joins two essential ligands, facilitating selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands linked together; one ligand targets an E3 ubiquitin ligase and the other targets the protein of interest. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6722 mL	18.3608 mL	36.7215 mL
5 mM	0.7344 mL	3.6722 mL	7.3443 mL
10 mM	0.3672 mL	1.8361 mL	3.6722 mL
50 mM	0.0734 mL	0.3672 mL	0.7344 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481