

Anti-FLT3 Antibody (AGS62P)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-FLT3 Antibody (AGS62P) is an ADC antibody targeting FLT3 for research in acute myeloid leukemia.
Targets(IC50)	FLT,ADC Antibody

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481