


PROTAC BRD4 ligand-3

Chemical Properties

CAS No. :	1445994-28-4
Formula:	C ₂₀ H ₁₅ F ₂ N ₃ O ₂
Molecular Weight:	367.35
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTAC BRD4 Ligand-3 serves as a ligand for target proteins in PROTACs (Ligands for Target Protein for PROTACs). It is utilized in the synthesis process.
Targets(IC50)	Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7222 mL	13.611 mL	27.222 mL
5 mM	0.5444 mL	2.7222 mL	5.4444 mL
10 mM	0.2722 mL	1.3611 mL	2.7222 mL
50 mM	0.0544 mL	0.2722 mL	0.5444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481