

PPO-IN-15

Chemical Properties


CAS No. : 3029518-71-3

Formula: C₁₇H₁₈ClFN₄O₅S

Molecular Weight: 444.87

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PPO-IN-15 (compound (R)-I-5) serves as an inhibitor of protoporphyrinogen IX oxidase (PPO), offering efficacy against resistant weeds while maintaining safety for wheat and rice.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2478 mL	11.2392 mL	22.4785 mL
5 mM	0.4496 mL	2.2478 mL	4.4957 mL
10 mM	0.2248 mL	1.1239 mL	2.2478 mL
50 mM	0.045 mL	0.2248 mL	0.4496 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481