

Anti-BRCA1 Antibody (BR64)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-BRCA1 Antibody (BR64) is a mouse IgG1, κ chimeric antibody that specifically targets human BRCA1. The recommended isotype control for the Anti-BRCA1 Antibody (BR64) is the Mouse IgG1 kappa, Isotype Control.
Targets(IC50)	DNA/RNA Synthesis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481