

Anti-NKG2D/CD314 Antibody (1D11)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-NKG2D/CD314 Antibody (1D11) is a chimeric mouse IgG1, κ targeting NKG2D/CD314 (human origin). Its isotype control is the MouseIgG1kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481