

Anti-Mouse Qa-1b Antibody (4C2.4A7.5H11)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-Mouse Nonclassical MHC Class I molecule Qa-1b Antibody (4C2.4A7.5H11) is a mouse-derived IgG1 antibody targeting the Qa-1b nonclassical MHC Class I molecule. The isotype control for this compound is Mouse IgG1 kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481