


Tyk2-IN-10

Chemical Properties

CAS No. :	2368945-72-4
Formula:	C ₂₅ H ₂₇ N ₅ O ₃
Molecular Weight:	445.51
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tyk2-IN-10 acts as an inhibitor of the Tyrosine Kinase 2 (Tyk2)-mediated signaling pathway, which plays a role in inflammation regulation.
Targets(IC50)	Others,JAK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2446 mL	11.2231 mL	22.4462 mL
5 mM	0.4489 mL	2.2446 mL	4.4892 mL
10 mM	0.2245 mL	1.1223 mL	2.2446 mL
50 mM	0.0449 mL	0.2245 mL	0.4489 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481