


## E3 Ligase Ligand-linker Conjugate 115

## Chemical Properties

CAS No. :	2375564-59-1
Formula:	C37H45FN4O7S
Molecular Weight:	708.84
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	E3 Ligase Ligand-linker Conjugate 115, a type of E3 ubiquitin ligase ligand-linker conjugate (E3 Ligase Ligand-Linker Conjugates), is utilized in the synthesis of PROTAC SMARCA2/4-degrader-27.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4108 mL	7.0538 mL	14.1076 mL
5 mM	0.2822 mL	1.4108 mL	2.8215 mL
10 mM	0.1411 mL	0.7054 mL	1.4108 mL
50 mM	0.0282 mL	0.1411 mL	0.2822 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481