


PROTAC SMARCA2/4-degrader-27

Chemical Properties

CAS No. :	2375564-54-6
Formula:	C49H58FN9O6S
Molecular Weight:	920.11
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PROTAC SMARCA2/4-degrader-27 (PROTAC 2) serves as a targeted degrader, utilizing PROTAC technology to degrade both SMARCA2 and SMARCA4.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0868 mL	5.4341 mL	10.8683 mL
5 mM	0.2174 mL	1.0868 mL	2.1737 mL
10 mM	0.1087 mL	0.5434 mL	1.0868 mL
50 mM	0.0217 mL	0.1087 mL	0.2174 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481