


## Sarubicin B

## Chemical Properties

CAS No. :	84745-01-7
Formula:	C <sub>13</sub> H <sub>10</sub> N <sub>2</sub> O <sub>4</sub>
Molecular Weight:	258.23
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sarubicin B, an active quinone antibiotic isolated from the culture filtrate of the Streptomyces strain (Streptomyces strain) JA 2861, effectively inhibits Gram-positive bacteria while showing ineffectiveness against Gram-negative organisms.
Targets(IC50)	Antibacterial, Antibiotic

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8725 mL	19.3626 mL	38.7252 mL
5 mM	0.7745 mL	3.8725 mL	7.745 mL
10 mM	0.3873 mL	1.9363 mL	3.8725 mL
50 mM	0.0775 mL	0.3873 mL	0.7745 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481