

E3 Ligase Ligand-linker Conjugate 117

Chemical Properties

CAS No. :

Formula: C34H44N8O3

Molecular Weight: 612.76

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3 Ligase Ligand-linker Conjugate 117, a conjugate of E3 ubiquitinase ligand and linker (E3 Ligase Ligand-Linker Conjugates), is utilized in the synthesis of PROTAC SMARCA2/4-degrader-29.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.632 mL	8.1598 mL	16.3196 mL
5 mM	0.3264 mL	1.632 mL	3.2639 mL
10 mM	0.1632 mL	0.816 mL	1.632 mL
50 mM	0.0326 mL	0.1632 mL	0.3264 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481