


NSC682769

Chemical Properties

CAS No. : 205048-79-9
 Formula: C₂₇H₂₂N₂O₃
 Molecular Weight: 422.48
 Storage: Keep away from direct sunlight
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	NSC682769, a PROTAC target protein ligand (Ligands for Target Protein for PROTACs), serves as a precursor in synthesis.
Targets(IC50)	Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.367 mL	11.8349 mL	23.6698 mL
5 mM	0.4734 mL	2.367 mL	4.734 mL
10 mM	0.2367 mL	1.1835 mL	2.367 mL
50 mM	0.0473 mL	0.2367 mL	0.4734 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481