


E3 ligase Ligand 33

Chemical Properties

| | |
|-------------------|--|
| CAS No. : | 3023410-43-4 |
| Formula: | C ₂₅ H ₃₄ N ₄ O ₄ |
| Molecular Weight: | 454.56 |
| Storage: | Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small> |



Biological Description

| | |
|---------------|---|
| Description | E3ligaseLigand 33 (56-9B) is an E3 ubiquitin ligase ligand utilized in cancer research. |
| Targets(IC50) | Ligands for E3 Ligase |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.1999 mL | 10.9996 mL | 21.9993 mL |
| 5 mM | 0.440 mL | 2.1999 mL | 4.3999 mL |
| 10 mM | 0.220 mL | 1.100 mL | 2.1999 mL |
| 50 mM | 0.044 mL | 0.220 mL | 0.440 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481