

## Anti-TROP-2 Antibody (Pr1E11)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	The Anti-TROP-2 Antibody (Pr1E11) is a chimeric mouse IgG1, κ antibody targeting TROP-2 (human). Its isotype control is Mouse IgG1 kappa, Isotype Control.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481