

Anti-Glycophorin A (type M) Antibody (6A7M)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Glycophorin A (type M) Antibody (6A7M) is a mouse IgG1, κ chimeric antibody designed to target human Glycophorin A (type M). The recommended isotype control for this antibody is Mouse IgG1 kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481