

E3 ligase Ligand 31

Chemical Properties

CAS No. : 2357114-13-5


Formula: C₁₆H₁₇N₃O₄

Molecular Weight: 315.32

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3ligaseLigand 31 is an E3 ubiquitin ligase ligand (Ligands for E3 Ligase). It is utilized in the synthesis of PTOTAC HSD17B13 degrader 1.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1714 mL	15.8569 mL	31.7138 mL
5 mM	0.6343 mL	3.1714 mL	6.3428 mL
10 mM	0.3171 mL	1.5857 mL	3.1714 mL
50 mM	0.0634 mL	0.3171 mL	0.6343 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481