


## SMARCA2/4-IN-2

## Chemical Properties

CAS No. :	1339378-93-6
Formula:	C <sub>10</sub> H <sub>9</sub> N <sub>3</sub> O
Molecular Weight:	187.2
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	SMARCA2/4-IN-2, serving as a PROTAC target protein ligand (Ligands for Target Protein for PROTACs), is utilized in synthesis.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3419 mL	26.7094 mL	53.4188 mL
5 mM	1.0684 mL	5.3419 mL	10.6838 mL
10 mM	0.5342 mL	2.6709 mL	5.3419 mL
50 mM	0.1068 mL	0.5342 mL	1.0684 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481