

## Anti-Fy3 Antibody (MIMA-29)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Fy3 Antibody (MIMA-29) is a mouse IgG2a, $\kappa$ chimeric antibody that targets human Fy3. The recommended isotype control for this antibody is Mouse IgG2a kappa, Isotype Control.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481