

HDAC-IN-77

Chemical Properties

CAS No. :

Formula: C₂₂H₂₆N₄O₂S

Molecular Weight: 410.53

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HDAC-IN-77 (HL-5s), an HDAC inhibitor, has the capability to induce ferroptosis and suppress the Nrf2/HO-1 signaling pathway. This compound is utilized in cancer research.
Targets(IC50)	Ferroptosis,Nrf2,HDAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4359 mL	12.1794 mL	24.3588 mL
5 mM	0.4872 mL	2.4359 mL	4.8718 mL
10 mM	0.2436 mL	1.2179 mL	2.4359 mL
50 mM	0.0487 mL	0.2436 mL	0.4872 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481