

Guanylin TFA (mouse, rat)

Chemical Properties

CAS No. :

Formula: C₆₀H₉₀N₁₆O₂₂S₄.xC₂H₅F₃O₂

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Guanylin (mouse, rat) TFA is a peptide composed of 15 amino acids and acts as an activator of intestinal guanylate cyclase. This compound is utilized in research related to diarrhea.
Targets(IC50)	Guanylate cyclase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481