

## Anti-IL-4 Antibody (MP4-25D2)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-IL-4 Antibody (MP4-25D2) is a rat IgG1, κ chimeric antibody that targets human IL-4. The recommended isotype control for Anti-IL-4 Antibody (MP4-25D2) is Rat IgG1 kappa, Isotype Control.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481