

## Anti-Rat FcRn heavy chain heterodimers Antibody (1G3)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-RatFcRnheavy chain heterodimers Antibody (1G3) is a mouse-derived antibody, specifically targeting the FcRn heavy chain heterodimers. The isotype control for Anti-RatFcRnheavy chain heterodimers Antibody (1G3) is designated as MouseIgG1kappa, Isotype Control.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481