

Anti-CD3 Antibody (OKT-3)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-CD3 Antibody (OKT-3) is a mouse IgG2a, κ chimeric antibody that targets human CD3. The recommended isotype control for Anti-CD3 Antibody (OKT-3) is the Mouse IgG2a kappa, Isotype Control.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481