

HPK1-IN-51

## Chemical Properties

CAS No. :

Formula: C<sub>28</sub>H<sub>31</sub>N<sub>5</sub>O

Molecular Weight: 453.58

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HPK1-IN-51, serving as a PROTAC target protein ligand (Ligands for Target Protein for PROTACs), is utilized in synthesis processes.
Targets(IC50)	Ligands for Target Protein for PROTAC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2047 mL	11.0234 mL	22.0468 mL
5 mM	0.4409 mL	2.2047 mL	4.4094 mL
10 mM	0.2205 mL	1.1023 mL	2.2047 mL
50 mM	0.0441 mL	0.2205 mL	0.4409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481