


## (S)-Dexfadrostat

## Chemical Properties

CAS No. :	102676-86-8
Formula:	C <sub>14</sub> H <sub>13</sub> N <sub>3</sub>
Molecular Weight:	223.27
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-Dexfadrostat ((S)-Fadrozole), an aromatase inhibitor, exhibits an IC <sub>50</sub> of 4.6 nM in in vitro assays using human placental microsomes. It is applicable for research on estrogen-dependent breast cancer, gynecomastia, and systemic lupus erythematosus.
Targets(IC <sub>50</sub> )	Others,Cytochromes P450

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4789 mL	22.3944 mL	44.7888 mL
5 mM	0.8958 mL	4.4789 mL	8.9578 mL
10 mM	0.4479 mL	2.2394 mL	4.4789 mL
50 mM	0.0896 mL	0.4479 mL	0.8958 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481