

Naphthyridin-Me-dimethoxybenzene-COOH

Chemical Properties

CAS No. :

Formula: C19H18N2O5

Molecular Weight: 354.36

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Naphthyridin-Me-dimethoxybenzene-COOH serves as a ligand for target proteins in PROTAC applications (Ligands for Target Protein for PROTACs). This compound is also utilized in synthetic processes.
Targets(IC50)	Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.822 mL	14.1099 mL	28.2199 mL
5 mM	0.5644 mL	2.822 mL	5.644 mL
10 mM	0.2822 mL	1.411 mL	2.822 mL
50 mM	0.0564 mL	0.2822 mL	0.5644 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481