

Pomalidomide-C2-NHCO-C6

Chemical Properties

CAS No. :

Formula: C₂₂H₂₈N₄O₅

Molecular Weight: 428.48

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pomalidomide-C2-NHCO-C6 serves as an E3 Ligase Ligand-Linker Conjugate. This compound facilitates the synthesis of PROTAC FGFR2 degrader 1.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3338 mL	11.6692 mL	23.3383 mL
5 mM	0.4668 mL	2.3338 mL	4.6677 mL
10 mM	0.2334 mL	1.1669 mL	2.3338 mL
50 mM	0.0467 mL	0.2334 mL	0.4668 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481