

Anti-SARS-CoV2 S protein Antibody(RBD epitope B,SARS2-34)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-SARS-CoV-2 S protein Antibody (RBD epitope B, SARS2-34) is a mouse-derived IgG1 kappa antibody that specifically targets the SARS-CoV-2 S protein. The recommended isotype control for this antibody is Mouse IgG1 kappa, Isotype Control.
Targets(IC50)	SARS-CoV

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481