


8-Hydroxyoctanoic acid-thalidomide

Chemical Properties

CAS No. :	2399454-65-8
Formula:	C ₂₁ H ₂₇ N ₃ O ₅
Molecular Weight:	401.46
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	8-Hydroxyoctanoic acid-thalidomide is a type of E3 ubiquitin ligase ligand-linker conjugate (E3 Ligase Ligand-Linker Conjugates). It is utilized in the synthesis of Lw13.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4909 mL	12.4545 mL	24.9091 mL
5 mM	0.4982 mL	2.4909 mL	4.9818 mL
10 mM	0.2491 mL	1.2455 mL	2.4909 mL
50 mM	0.0498 mL	0.2491 mL	0.4982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481