

## E3 Ligase Ligand-linker Conjugate 116

## Chemical Properties

CAS No. : 2409844-90-0


Formula: C<sub>48</sub>H<sub>75</sub>N<sub>5</sub>O<sub>15</sub>S

Molecular Weight: 994.2

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	E3 Ligase Ligand-linker Conjugate 116, an E3 ubiquitinase ligand plus linker conjugate (E3 Ligase Ligand-Linker Conjugates), facilitates the synthesis of PROTAC SMARCA2/4-degrader-28.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0058 mL	5.0292 mL	10.0583 mL
5 mM	0.2012 mL	1.0058 mL	2.0117 mL
10 mM	0.1006 mL	0.5029 mL	1.0058 mL
50 mM	0.0201 mL	0.1006 mL	0.2012 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481