

Anti-Mouse TIM-3 Antibody (RMT3-23)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-MouseTIM-3Antibody (RMT3-23) is a rat-derived IgG2akappa antibody targeting TIM-3 in mice. The isotype control for Anti-MouseTIM-3Antibody (RMT3-23) is RatIgG2akappa, Isotype Control.
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481