

(R,S,R)-AHPC-PEG2-C2-boc

Chemical Properties

CAS No. :

Formula: C35H52N4O8S

Molecular Weight: 688.87

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(R,S,R)-AHPC-PEG2-C2-boc is classified as an E3 Ligase Ligand-Linker Conjugate. This compound is instrumental in the synthesis of PROTAC YAP degrader-1.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4517 mL	7.2583 mL	14.5165 mL
5 mM	0.2903 mL	1.4517 mL	2.9033 mL
10 mM	0.1452 mL	0.7258 mL	1.4517 mL
50 mM	0.029 mL	0.1452 mL	0.2903 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481