


## PROTAC SMARCA2/4-degrader-28

## Chemical Properties

CAS No. :	2409844-88-6
Formula:	C <sub>54</sub> H <sub>68</sub> ClN <sub>9</sub> O <sub>11</sub> S <sub>2</sub>
Molecular Weight:	1118.75
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



## Biological Description

Description	PROTAC SMARCA2/4-degrader-28 (PROTAC 1) functions as a partial degrader of SMARCA2 and SMARCA4 through the PROTAC-based mechanism.
Targets(IC50)	Epigenetic Reader Domain,PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8939 mL	4.4693 mL	8.9385 mL
5 mM	0.1788 mL	0.8939 mL	1.7877 mL
10 mM	0.0894 mL	0.4469 mL	0.8939 mL
50 mM	0.0179 mL	0.0894 mL	0.1788 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481