


Androgen receptor degrader-5

Chemical Properties

CAS No. :	2902679-11-0
Formula:	C ₂₉ H ₂₅ F ₄ N ₅ O ₂
Molecular Weight:	551.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Androgen receptor degrader-5 (compound 14k) demonstrates superior capabilities, specifically in androgen receptor degradation and antiproliferative activity, showcasing its promising properties.
Targets(IC50)	Androgen Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8131 mL	9.0657 mL	18.1314 mL
5 mM	0.3626 mL	1.8131 mL	3.6263 mL
10 mM	0.1813 mL	0.9066 mL	1.8131 mL
50 mM	0.0363 mL	0.1813 mL	0.3626 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481