

PKMYT1-IN-3

Chemical Properties

CAS No. :

Formula: C₂₄H₂₆FN₅O₂

Molecular Weight: 435.49

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PKMYT1-IN-3 (compound 8ma) is a potent, selective inhibitor of PKMYT1, demonstrating an IC ₅₀ of 16.5 nM and exhibiting antitumor activity.
Targets(IC ₅₀)	Wee1

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2963 mL	11.4813 mL	22.9626 mL
5 mM	0.4593 mL	2.2963 mL	4.5925 mL
10 mM	0.2296 mL	1.1481 mL	2.2963 mL
50 mM	0.0459 mL	0.2296 mL	0.4593 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481