

Anti-LFA-1 α /CD11a Antibody (TS-1/22.1.1.13)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-LFA-1 α /CD11a Antibody (TS-1/22.1.1.13) is a chimeric mouse IgG1 antibody targeting LFA-1 α /CD11a of human origin. Its isotype control is designated as MouseIgG1kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481