

## Anti-MUC1/CD227 Antibody (C595 (NCRC48))

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	The Anti-MUC1/CD227 Antibody (C595 (NCRC48)) is a chimeric mouse IgG3, κ antibody targeting human MUC1/CD227. Its isotype control is the Mouse IgG3 kappa, Isotype Control.
Targets(IC50)	Mucin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481