

## Anti-Mouse IL-4 Antibody (BVD6-24G2)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-MouseIL-4Antibody (BVD6-24G2) is a rat-derived IgG1kappa antibody that targets IL-4 in mice. The isotype control for Anti-MouseIL-4Antibody (BVD6-24G2) is RatIgG1kappa, Isotype Control.
Targets(IC50)	Interleukin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481