

ABBV-383

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

|               |                                                                                                                                                                                                                                                                                                                                                  |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description   | ABBV-383 (ABBV-383), a BCMA × CD3 bispecific T-cell engager (BiTE), targets and inhibits the B-cell maturation antigen (BCMA) while concurrently activating the T-cell surface glycoprotein CD3 complex. This compound is utilized in research focused on multiple myeloma, immunoglobulin light chain amyloidosis, and cardiovascular diseases. |
| Targets(IC50) | Integrin                                                                                                                                                                                                                                                                                                                                         |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481