

Efpegarglucagon

Chemical Properties

CAS No. : 2766774-81-4

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Efpegarglucagon comprises the Fc fragment of human immunoglobulin G4 (IgG4), forming a dimer linked by disulfide bridges. One of the Fc chains is connected to a glucagon-like peptide through a linker (maleimide-polyethylene glycol-aldehyde).
Targets(IC50)	Glucagon Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481