

Anti-GPC3 Antibody (YP7)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-GPC3 Antibody (YP7) is a chimeric mouse IgG1, κ antibody that specifically targets human GPC3. The recommended isotype control for Anti-GPC3 Antibody (YP7) is Mouse IgG1 kappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481