


GSI-18

Chemical Properties

CAS No. : 362653-62-1
 Formula: C₁₇H₁₉NO₂S₂
 Molecular Weight: 333.47
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	GSI-18, a γ -secretase inhibitor, exhibits anticancer activity by mitigating anchorage-independent growth in pancreatic cancer cells through the inhibition of Notch signaling.
Targets(IC50)	Gamma-secretase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9988 mL	14.9939 mL	29.9877 mL
5 mM	0.5998 mL	2.9988 mL	5.9975 mL
10 mM	0.2999 mL	1.4994 mL	2.9988 mL
50 mM	0.060 mL	0.2999 mL	0.5998 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481