


2-Epitripdiolide

Chemical Properties

CAS No. : 74409-90-8
 Formula: C₂₀H₂₄O₇
 Molecular Weight: 376.4
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Epitripdiolide, a diterpenoid compound derived from <i>Tripterygium hypoglaucum</i> , exhibits notable antitumor activity. The IC ₅₀ values for A549, KB, KBvin, and MDA-MB-231 cell lines are 1.34, 0.4, 4.4, and 1.2 nM, respectively.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6567 mL	13.2837 mL	26.5675 mL
5 mM	0.5313 mL	2.6567 mL	5.3135 mL
10 mM	0.2657 mL	1.3284 mL	2.6567 mL
50 mM	0.0531 mL	0.2657 mL	0.5313 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481