

C1q Protein (Mouse)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	C1qProtein (Mouse) is a subunit of the complement system. It recognizes a variety of ligands and plays roles in maintaining immune system balance, promoting tissue repair, and defending against pathogens.
Targets(IC50)	Complement System

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481