

Anti-LCMV Nucleoprotein Antibody (VL-4)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	The Anti-LCMV Nucleoprotein Antibody (VL-4) is a mouse-derived IgG2akappa antibody that targets the LCMV nucleoprotein and originates from rats. The isotype control for this antibody is RatIgG2akappa, Isotype Control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481